

DNA Banking

Referred Out Test:

DNA Retrieval for Referred Out Test

C340

(*Note: Related charges may be billed back to the institution)		
SPECIMEN REQUIREMENTS		
 Fetal Tissue for Genetic Testing: Product of Conception (POC): 25 mg of UNFIXED tissue in sterile transport media or saline. Umbilical cord: 1 cm piece of UNFIXED tissue in sterile transport media or saline Skin / Cartilage: 0.5 cm x 0.5 cm piece of UNFIXED tissue in sterile transport media or saline 	 Paraffin Embedded Tissue for Aneuploidy Detection/Genotyping: Placental 6 plain slides sequentially cut of 10 microns thickness plus one H&E slide with fetal and maternal areas clearly marked Fetal 6 plain slides sequentially cut of 10 microns. 	
INSTRUCTIONS FOR SUBMISSION OF SPECIMENS		
Deliver specimens to the Laboratory as soon as possible. Maintain samples at room temperature during transit (do not refrigerate or freeze).		

MRN: _____

Health Card #:_____

Date of Birth:

For more details on collection, shipping and reference information , click here to access our Test Catalogue: https://www.sinaihealth.ca/areas-of-care/pathology-and-laboratory-medicine